



An Intelligent, Endoscopic Solution for Obesity

Introducing the first and only
LONG-TERM & REVERSIBLE
endoscopic obesity solution, with the
power to address
obese and super-obese patients.

ASPIRE

The Power of the AspireAssist®




[34]
kgs lost

[28]
kgs lost

[63]
kgs lost

[65]
kgs lost





These four siblings/in-laws lost 190 kgs collectively in 10 months

AspireAssist: Dual Mechanism for Long-Term Results

30 Percent
CALORIC REDUCTION

Removes a portion of
stomach contents,
bypassing absorption



Integrated
INTAKE MODULATION

6mm tube diameter re-trains
patient to chew extensively,
naturally slowing intake

Primary Mechanism: Caloric Reduction



Allow 20 minutes
after each meal for
digestion



In privacy of restroom,
remove AspireAssist®
from pocket-sized bag



Removes up to 30%
of calories, directly
into toilet



Discreet port "button"
when not in use


Patients aspirate after each main meal, three times per day.

Secondary Mechanism: Intake Modulation



Increased Chewing


(May/June 2015)



I chew my food more thoroughly since beginning therapy with the AspireAssist.


Increased Meal Consumption Time

(May/June 2015)



Reduced Calorie Consumption

(May/June 2015)




Since beginning therapy with the AspireAssist, I believe the number of daily calories I eat has:

Consumption is modulated by the 6mm tube diameter.

Patients must chew thoroughly to facilitate aspiration, training patient to eat more slowly for long-term weight maintenance.


Lifestyle Counseling Reinforces Healthy Habits



All patients receive lifestyle counseling with AspireAssist therapy, typically about 10 sessions in first year


Proven Clinical Results

Mean 2-Year Weight Loss



Comparable Results to Bariatric Surgery...

Mean % Excess Weight Loss (One Year)



...with the safety profile for even high-BMI patients

Excellent Track Record of Safety



NO SERIOUS COMPLICATIONS AFTER >300 PATIENTS

DRAWS ON 35 YEARS OF PEG TUBE BEST PRACTICES

EXCELLENT PATIENT TOLERANCE

NO ELECTROLYTE OR METABOLIC IMBALANCES

Safe and Routine Procedure

STANDARD PEG PROCEDURE REQUIRES NO ADDITIONAL TRAINING

The AspireAssist is implanted using the Ponsky “pull” PEG technique, a simple and routine procedure for Gastroenterologists & Surgeons

PROCEDURE FEATURES

15-MINUTE PROCEDURE

to place the tube endoscopically through the mouth

USES CONCIIOUS SEDATION


General anesthesia is typically not necessary

OUTPATIENT

Patients return home within 2 hours

HIGH PROCEDURE SUCCESS

Adequate transillumination in >99% of patients



Long-term and Reversible. Best of Both Worlds.

The AspireAssist is **the only endoscopic solution** that is both long-term and reversible, putting the patient in the driver's seat.

LONG-TERM SOLUTION

Intended for long-term use, although patients typically reduce frequency of use as they approach goal weight and adopt healthier habits


EASILY REVERSIBLE

Tube can be removed at any time in a simple 10-minute endoscopy

Patient Acceptability

Strong Patient Interest in AspireAssist



(International Market Research 2014)



Based on the information provided, please rate your overall level of interest in Therapy X.

Over 90% Satisfaction Rate Among US Study Participants

(US Pivotal Trial Survey, 2015)



Multiple European Studies Ongoing

Ongoing studies across Europe continue to support excellent weight loss results and safety profile



Swedish Post-Market Study

- ❖ 25 Subjects
- ❖ BMI 35-49

6-Month Results Published in *Endoscopy* 2015:
41% EWL, 36 pounds

Pan-European Post-Market Registry

- ❖ Enrolling 50 Subjects
- ❖ Spain, Austria, Czech, Italy, Belgium, Greece, UK
- ❖ BMI 35-65

Germany & Austria Post-Market Study

- ❖ Enrolling 30 Subjects
- ❖ BMI 35-65

Head-to-Head vs Gastric Bypass

- ❖ Enrolling 100 Subjects
- ❖ Sweden
- ❖ BMI 35-49

Super-Obese

- ❖ Enrolling 30 Subjects
- ❖ Czech Republic, Spain, Belgium, France
- ❖ BMI 59-79 to date

PATHWAY U.S. Pivotal Trial Design

171 SUBJECT TRIAL ACROSS 10 LEADING INSTITUTIONS

- Body Mass Index (BMI) 35 – 55
- Failed previous weight loss attempts

2:1 RANDOMIZATION

- 111 AspireAssist, 60 Lifestyle Therapy

PRIMARY ENDPOINTS

- Mean percent Excess Weight Loss (EWL) >10% over control at 52 weeks
- At least 50% “Responder Rate” at 52-weeks (defined as 25% EWL)


Institutions

Boston Medical Center
Brigham & Women’s Hospital
Weill Cornell Medical College
St. Mary Medical Center
University of Pennsylvania
Howard University
Northwestern University
Mayo Clinic
Washington University
VA Center/ UC San Diego

Data locked in June 2015

BOTH PRIMARY ENDPOINTS WERE MET; TRIAL SUCCESSFUL

Entering Final Stage of FDA Approval Process



Expected U.S. launch in early-to-mid 2016

ASPire
assist®