



LATITUDE™ AIR-CHARGED SINGLE-USE CATHETER

Simpler Than Solid State

Air-Charged Catheters make motility testing simpler than ever before. Disposable Latitude™ catheters set up in seconds and protect against cross-contamination. Latitude catheters are safe, easy to use, and more comfortable for patients. Latitude frees you from many encumbrances that complicate your procedures, slow you down, and keep you from expanding your practice. For advanced performance made easy, the choice is simple.

Air-Charged Catheters: Simpler than Solid State:

- Disposable catheter: Clean, single-use catheters.
- No cross-contamination risk: Every patient gets a new catheter.
- No repair or sterilisation costs: No downtime or backup catheters.
- Enhanced patient tolerance: Small catheter maximises comfort.
- Easier insertion articulating tip aids in difficult insertions.
- Every sensor circumferential: Eliminates anatomical irregularities.
- Low capital cost: No expensive capital outlay.



The Latitude Calibrator is available with easy-to-read gauge.



Pre-attached balloon makes for swift setup. Small diameter eliminates waiting for sphincter acclimation.



Optional, reusable stylet gives you a choice of insertion stiffness!



Advanced technology is simple to operate. Colour coded luer connection ensure proper setup.



Articulating tip aids in difficult insertions and the small catheter diameter maximises patient comfort, procedure duration and results.

SynMed Ltd.
SynMed House
7 The Pavilion Business Centre
6 Kinetic Crescent
Innova Park
Enfield
EN3 7FJ
Tel: +(44) 01992 782570
Fax: +(44) 01992 667010
www.synmed.co.uk
sales@synmed.co.uk

Manufactured by:
Clinical Innovations, Inc.
www.clinicalinnovations.com

CE
0344